

ICBNh 5173

Integrated fridge-freezer with
BioFresh Professional and NoFrost
Δ Peak

LIEBHERR

Fully Integrated



Fully Integrated

ICBNh 5173

Integrated fridge-freezer with
BioFresh Professional and NoFrost
△ Peak

↑
178
↓



LIEBHERR

G Made
in Germany



Bio Fresh	Bio Fresh Professional	Hydro Breeze	No Frost	Duo Cooling	Light Tower	Touch & Swipe-Display	Glass & Steel Interior	Ice Maker	Soft System
Soft-Telescopic Rails	Vario Space	Super Silent							

Type	Fixed door / Built-in appliance for integrated use
Total volume	246 l ¹
Total net	Fridge: 172 l / Freezer compartment: 67 l
Noise output	34 dB(A)
Frequency / voltage	50-60 Hz / 220-240 V~
Connection Rating	1.4 A 183 W
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55
Controls	
Type of Control	2.4" TFT-colour display, touch- & swipe Digital temperature display for fridge and freezer compartments
Interface	SmartDeviceBox: optional accessory (Only available in selected countries)
Door alarm	visual and acoustic
Alarm in case of malfunction	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with Ambilight and FreshAir activated charcoal filter
Door interior	GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray
Bottle shelf	Integrated bottle shelf
Interior light	LED light panel, both sides LED BioFresh lighting
BioFresh-Safes	2 BioFresh Safes with SoftTelescopic, 1 FlexSystem
Freezer compartment	
Star coding	****
Interior light	LED lighting
Number of drawers	3 drawers, of which 2 with telescopic partial pull-out, VarioSpace
IceMaker	with fixed 3/4" water connection, LED lighting
Ice cube capacity in 24 hours	1.2 kg
Ice cube stock	5.0 kg
Freezing capacity in 24 hours	6 kg
Features	
Door hinges	Reversible door hinge, left
Replaceable door seal	yes
EAN-Nr.	4016803115274

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

Features



BioFresh Professional

The BioFresh Professional safe gives food an additional boost of freshness for a long shelf life – wow effect included: Thanks to HydroBreeze a cold mist covers fruit and vegetables like a protective balm. In combination with a temperature just above 0 °C and high humidity, unpackaged fruits and vegetables feel right at home.



BioFresh with HydroBreeze

Do you want to cool fruit and vegetables like professionals? You'll be inspired by HydroBreeze. The cold mist combined with a temperature in the safe of just over 0 °C gives food that extra boost for longer shelf life. And it creates a genuine visual wow effect. HydroBreeze is activated every 90 minutes for 4 seconds and for 8 seconds when the door is opened.



Meat & Dairy safe

Do you want the tender veal fillet or fresh milk to stay fresh for as long as possible? The BioFresh safes will ensure this is the case. In the Meat & Dairy Safe, the temperature is close to 0 °C – ideal storage conditions for even the most delicate of foods. And best of all, everything is perfectly adjusted in the factory for immediate use – and for longer freshness.



NoFrost

When you open your freezer, you want to see frozen food – and certainly not ice and frost. NoFrost protects the freezer compartment from unwanted icing, which consumes a lot of energy and can be expensive. NoFrost means no more tedious and time-consuming defrosting of the freezer compartment, more time for other things – and saving money.



DuoCooling


DuoCooling uses two completely separate refrigeration circuits to ensure no air is exchanged between the fridge and freezer compartments. The food neither dries out nor is odour transferred. This means throwing less away, shopping less often, but saving more and greater enjoyment.




Ambient lighting

Our ambient lighting brings a special design touch to your refrigerator. The all-around lighting produces unique effects in the interior, enhancing the overall atmosphere. Paired with the SmartSteel-designed rear wall and the sleek black high-gloss fan, it creates a captivating and ambient space.

Technical drawing



Maximum weight of furniture door:
Fridge compartment: 18 kg
Freezer compartment: 12 kg



Products and their details are staged for advertising purposes. Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any damage. (20.03.2025)



home.liebherr.com