

Sender	Receiver							
	Germany Austria Colmar, France		Switzerland		Europe		Liebherr companies worldwide	
Shipment weight	≤100 kg*	>100 kg	≤100 kg*	>100 kg	≤100 kg*	>100 kg	≤100 kg*	>100 kg
Belgium	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Bulgaria	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Denmark	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Germany	TNT	Fried-Sped Noerpel	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Estonia	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Finland	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
France	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Greece	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Ireland	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Italy	TNT	Fried-Sped Gebrüder Weiss (only EMT & LLG)	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Croatia	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Latvia	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Lithuania	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Luxembourg	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Netherlands	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped

*Shipment weight up to 100kg or maximum package dimensions of 2.4m x 1.2m x 1.8m

** Please use the shipping instruction "Routing-Order Switzerland" for all collections in Switzerland

Sender	Receiver							
	Germany Austria Colmar, France		Switzerland		Europe		Liebherr companies worldwide	
Shipment weight	≤100 kg*	>100 kg	≤100 kg*	>100 kg	≤100 kg*	>100 kg	≤100 kg*	>100 kg
Norway	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Austria	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Poland	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Portugal	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Romania	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Sweden	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Switzerland**	Fried-Sped	Fried-Sped	≤30 kg Schweizerische Post	>30 kg Zumwald	Fried-Sped	Fried-Sped	Fried-Sped	Fried-Sped
Slovakia	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Slovenia	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Spain	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Czech Republic	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Hungary	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
United Kingdom	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped	TNT	Fried-Sped
Countries not listed	Please contact your responsible Liebherr company							

*Shipment weight up to 100kg or maximum package dimensions of 2.4m x 1.2m x 1.8m

** Please use the shipping instruction "Routing-Order Switzerland" for all collections in Switzerland