

EFI 2853


Impulse Sales Chest Freezer



LIEBHERR

EFI 2853

Impulse Sales Chest Freezer




home.liebherr.com



EFI 2853	
Total gross/net capacity l (cu.ft.)	9.7 / 7.0 (276 / 197)
Exterior dimensions in inches (w/d/h) in mm (w/d/h)	41 3/16" / 26 13/16" / 32 1/2" 1045 / 680 / 825
Interior dimensions in inches (w/d/h) in mm (w/d/h)	36 1/4" / 21 1/8" / 19 5/16" 920 / 535 / 490
Total display area in m² (TDA)/ft²	0.54 / 5.81
Energy Consumption in 365 days*	593 kWh
Ambient temperature in °C (°F)	50°F to 86°F (10°C to 30°C)
Refrigerant	R290
Sound rating**	42 dB(A)
Connection rating	2.5 A/219 W
Frequency/voltage	60 Hz/115 - 127V~
Cooling system/Defrost	Static/Manual
Temperature range in °C (°F)	14°F to -11°F (-10°C to -24°C)
Lid material	Sliding glass lid
External cabinet finish	Steel/White
Internal liner material	Steel
Frame color	Grey
Type of control / Temperature display	Analog control/internal analogue
Interior Light	LED
Number of baskets	4
Number of baskets, max. (see picture)	4
Castors	Swivel casters
Lock	Can be retrofitted
Insulation in inches (mm)	2 3/8" (60)
Water drain	Yes
Gross/net weight in Kg (lbs)	135.1 / 118.6 (61.3 / 53.8)

Unit Dimensions



Accessories

Basket 10 5/8 in (270 mm)	7111009
Basket 8 1/4 in (210 mm)	7111039
Partition grid 10 5/8 in (270 mm) / High grid	7111031
Partition grid 10 5/8 in (270 mm) / Low grid	7111033
Partition grid 8 1/4 in (210 mm)	7111787
Rubber bumpers, gray	7266707
Cylinder lock for glass sliding lid	7043359

* Measured unloaded at 26.6°C (80°F) ambient temperature with set point of 3.3°C (38°F).

** Noise level – sound power re 1pW – according to IEC 60704-2-14.

LIEBHERR

Subject to modification. For up-to-date information, visit home.liebherr.com.