


Superb Conditions  
for Fine Wines

Quality, Design and Innovation



Wkb 1712

**LIEBHERR**



**Optimum humidity:** Liebherr wine cabinets guarantee an optimum humidity and this ensured the corks are supple in place whilst still allowing the wine to breathe in its maturing process.



**LED lighting:** The tinted glass door provides not only a perfect UV protection, but also an optimum glance at wine. The integrated LED lighting ensures a pleasant illumination of the complete interior.




**Wooden shelves:** The wooden shelves are ideal for storing Bordeaux bottles. They can be stored in opposite directions to make full use of the capacity.

The **WKb 1712** wine cabinet, our most compact model from the Liebherr Barrique series, can store up to 60 Bordeaux bottles. This unit ensures a constant temperature that can be set between 5°C to 20°C. Wines are stored optimally on three untreated height-adjustable wooden shelves, all protected from UV radiation thanks to a tinted insulated glass door. The LED lighting ensures that the interior is illuminated evenly and pleasantly so wines can be presented with illumination even over a longer period.

Optional accessories<sup>1</sup> include our charcoal filter that protects your finest collection against odour; best changed once every 12 months. For added flexibility, opt for our add-on shelves.

**WKb 1712** Barrique LED Wine

Sound rating dB(A) *	40
Total capacity litre	121
Max No. of Bordeaux bottles (0.75l) **	60
Number of shelves	3
Shelf material	Beech wood
Electronic control	LCD display, buttons
Temperature range °C	+5°C to +20°C
Temperature zones	1
Temperature alarm	Visual and audible
Interior light	LED
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Tinted glazed door, black frame
Integrated lock	Yes
SoftSystem	No
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	85 / 65.7 / 66.3
Gross / Net weight kg	50 / 45
Voltage V/Hz	220-240 / 50



Activated Charcoal Filter  
Part no.: 988129100

\* Noise output - according to EN 60704-3

\*\* These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

<sup>1</sup> Available accessories: Activated Charcoal Filter (Part no. 988129100), Regular Shelf (Part no. 711173400)

Note: Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any loss.