

EFE 4652




LIEBHERR


Impulse Sales Chest Freezer



| Impulse Sales Chest Freezer | EFE 4652 | LED |
|---|---|-----|
| Total gross/net capacity cu.ft. (l) | 15.7 / 12.0 (444 / 340) | |
| Exterior dimensions in inches (w/d/h) in mm (w/d/h) | 57 11/16" / 26 13/16" / 32 1/2" 1465 / 680 / 825 | |
| Interior dimensions in inches (w/d/h) in mm (w/d/h) | 52 13/16" / 21 1/8" / 20 5/8" 1340 / 535 / 523 | |
| Total display area in m² (TDA) / ft² | 0.78 / 8.39 | |
| Energy consumption in 365 days * | 811 kWh | |
| Ambient temperature in °F (°C) | 50°F to 95°F (10°C to 35°C) | |
| Refrigerant | R290 | |
| Sound rating ** | 52 dB(A) | |
| Connection rating | 2.5 A/238 W | |
| Frequency / Voltage | 60 Hz/115 – 127V~ | |
| Cooling system | static | |
| Defrost | manual | |
| Temperature range °F (°C) | 14°F to -11°F (-10°C to -24°C) | |
| Lid material | sliding glass lid | |
| External cabinet finish | steel/white | |
| Internal liner material | steel | |
| Frame color | grey | |
| Type of control | analog control | |
| Temperature display | internal analogue | |
| Interior light | LED | |
| Number of baskets | 6 | |
| Number of baskets, max (see picture) | 6 | |
| Castors | swivel casters | |
| Lock | can be retrofitted | |
| Insulation in inches (mm) | 2 3/8" (60) | |
| Water drain | yes | |
| Gross/net weight in lbs (kg) | 172 / 152.3 (78 / 69.1) | |
| Accessories | | |
| Basket 270 mm | 7111095 | |
| Basket 210 mm | 7111097 | |
| Partition grid 270 mm | 7111029 | |
| Partition grid 210 mm | 7111784 | |
| Rubber bumpers, grey | 7266707 | |
| Cylinder lock for glass sliding lid | 7043660 | |
| Scoop ice cream containers | on request | |
| Baskets for scoop ice cream containers | on request | |
| External thermostat dial | on request | |



Unit Dimensions



Subject to modification. For up-to-date information, visit home.liebherr.com.