



ENERGY STAR CERTIFIED

# Lab Grade Refrigerators and Freezers

## Liebherr - LRT50W2HC : LRT50W2HC

### Specifications

ENERGY STAR Unique ID:	2323102
ENERGY STAR Partner:	Liebherr-Hausgerate Ochsenhausen GmbH
Brand Name:	Liebherr
Model Name:	LRT50W2HC
Model Number:	LRT50W2HC
Product Type:	High Performance Refrigerator
Total Volume (cubic feet):	46.79
Lab Grade Refrigerator or Freezer Energy Consumption (kWh/day):	3.65
Peak Temperature Variance (°C):	2.3
Product Form Factor:	Upright
Total Number of Doors:	2
Number of glass doors:	0
Number of solid doors:	2
Door Opening Orientation:	Swinging
Number of Drawers or Shelves:	8
Number of Internal Ports and/or Access Holes:	1
Depth (in):	32.88
Height:	85.06
Width:	56.32
Defrost Type:	Automatic
Refrigerant with GWP:	R-290 (GWP:3   Lower GWP)
Installed Accessories:	None
Date Available on Market:	1514782800000
Date Certified:	1531886400000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Captured On:  
06/10/2024