

# Bleed Test Set 99127B03

## For Engine Bleed Air System Trouble Shooting



### • BTS Users

Liebherr-Aerospace is proud to deliver its 25th Bleed Test Set (BTS) to the growing community of BTS users amongst the Airbus operators.

The BTS is a compact and mobile digest means. It allows to decrease significantly the unjustified component removals rates and enhances the trouble shooting of bleed systems and components on the Airbus Family

aircraft from A318 to A380. It can be a key tool in enhancing aircraft bleed system operational reliability.

Officially recognized and endorsed by Airbus in the aircraft maintenance manual for the single aisle and long range fleet, the BTS is a special test equipment which can be used for on-wing or off-wing maintenance at various stages of the component's life.

# LIEBHERR

# The benefits of using the Bleed Test Set are:

- Improved trouble shooting accuracy and time
- Testing of pneumatic components during trouble shooting and isolation of individual component failures
- Health check program application using OEM or AMM documentation for preventative maintenance
- Reduction of unjustified removals and unnecessary returns to maintenance provider workshops
- Access to a bleed sense line tester included in the Bleed Test Set



ATA36 Engine Bleed Air System

