




ELEVATION LEFT WTes 5972
1 : 10



ELEVATION FRONT WTes 5972
1 : 10



FLOOR PLAN WTes 5972
1 : 5



3D View WTes 5972

WTes 5972 Vinidor	
Total capacity litre	526
Capacity of wine compartment litre	526
Max. No. of Bordeaux btls 0.75	211 *
Shipping weight kg	137
Climate rating	T
Door hinges	right, reversible
Outside dimensions in cm (H / W / D)	192 / 70 / 74.2
Energy supply Voltage V/Hz	220 – 240 / 50
Electronic control	LCD display in the interior, touch electronic Digital temperature display for both wine zones

* These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

WTes 5972_136

Scale

As indicated