




ELEVATION LEFT WTes 5872
1 : 10



ELEVATION FRONT WTes 5872
1 : 10



FLOOR PLAN WTes 5872
1 : 5



3D View WTes 5872

WTes 5872 Vinidor	
Total capacity litre	511
Capacity of wine compartment litre	511
Max. No. of Bordeaux btls 0.75	178 *
Shipping weight kg	148
Climate rating	T
Door hinges	right, reversible
Outside dimensions in cm (H / W / D)	192 / 70 / 74.2
Energy supply Voltage V/Hz	220 – 240 / 50
Electronic control	LCD display in the interior, touch electronic Digital temperature display for 3 wine zones

* These models are only bound for storing wine. Maximum DIN norm
bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm,
Ø 76.1 mm)

WTes 5872_136

Scale

As indicated