

---

# Equipment for concrete pumps

---

Equipment catalogue for safe and  
comfortable working

## LIEBHERR

Concrete pumps



# Our equipment make the difference

## General information

Here you will find an overview of the equipment and accessories of our Liebherr truck-mounted concrete pumps. Our range of additional equipment covers a wide range of applications.

Customers benefit from a wide range of additional accessories from all areas and customisation.

**The illustrations are not binding. There are variations depending on the model and design.**



You can also find all our  
accessories online at:  
[www.my.liebherr.de](http://www.my.liebherr.de)

QR code MyLiebherr:



# Cleaning


## Drip bag

Item number: 14301180




## Coupling adapter for water tank connection 52-C to 3/4" GHT

Item number: 14933058



## Coupling adapter for water tank connection C - GEKA


Item number: 14273234






## Coupling adapter for water tank connection C - CAMLOCK

Item number: 14273170



## Coupling adapter for water tank connection C - DSP

Item number: 14273152



## Low-pressure cleaning gun

Item number: 13688860



## PowerControl jet pipe stepless 800 mm

Pressure adjustment via rotary handle

Manufacturer: Kärcher

Item number: 13911731

## Cleaning pipe 400 mm

For information: cleaning pipe 200 mm standard version

Manufacturer: Kärcher

Item number: 13911729


## HD rotor nozzle (dirt cutter)

Manufacturer: Kärcher

Item number: 14399932



Please use only in  
conjunction with  
matched Kärcher  
components.



# Storage options truck mounted concrete pumps


**!** Please contact before installation your responsible service partner.

## Retrofit kit storage tray aluminum (assembly on left leg)

**Length:** 4 m  
**Capacity:** 20 m

**Ø hoses:** 125 mm  
**Material:** Aluminium

**Item number:** 13854097



## Retrofit kit hose & pipe tray (to be mounted on left inside leg)

**Tubes:** 6 pieces  
**Ø Tubes:** 100 mm  
**Material:** PE plastic

**Length:** 3 m  
**Capacity:** 18 m

**Item number:** 14272798



## Retrofit kit hose holder with railing (to be mounted on right leg)

**Hoses:** 6 pieces  
**Ø Hoses:** 100 mm  
**Material:** PE plastic

**Length:** 5 m  
**Capacity:** 30 m

**Item number:** 14395069



## Retrofit kit hose holder (to be mounted on left leg)

Hoses: 6 pieces

Ø Hoses: 100 mm

Material: PE plastic

Length: 5 m

Capacity: 30 m


Item number: 14289606



## Stainless steel tool box

Dimensions: 800x500x500 mm


Item number: 13598995



## Toolbox PE

Dimensions: 800x400x400 mm

Item number: 13752536






## Toolbox PE

Dimensions: 500x700x60 mm


Item number: 12681313



## Toolbox PE

Dimensions: 388x650x600 mm


Item number: 12679663



## Closed document container

Dimensions: 250x100x450 mm


Item number: 13787239



# Pipe container

Dimensions: 205x205x600 mm

Item number: 14429218




# Storage options crawler concrete pumps

**!** Please contact  
before installation  
your responsible  
service partner.

## Tool box incl. holder PE

Dimensions: 800x400x400 mm


Item number: 14447015



## 1 x toolbox (PE) incl. holder

Dimensions: 400x370x350 mm


Item number: 14429389



## 2 x toolbox (PE) incl. holder

Dimensions: 2x 400x370x350 mm

Item number: 14429384



## Toolbox (PE)

**Dimensions:** 500x400x370 mm


**Item number:** 14442265



## Closed document container

**Dimensions:** 250x100x450 mm

**Item number:** 14448125





# Radio remote control

## Y braces supplemented to the hip support set

Item number: 13681420



## Cross shoulder strap

Item number: 13681422



## Shoulder hook, technos B

Item number: 11635235



## Hip belt for technos B

Item number: 10181431



## Hip belt with adjustable transmitter mounting technos B

Item number: 13681429



## Universal hip belt with sliding transmitter mounting

Item number: 13681419





## Extension for waist belt

Item number: 12697567



## Protective bag for technos B

Item number: 13681418



## Storage device for technos B (to be mounted inside truck cabine)

Item number: 12275333



## **Storage device for technos B (to be mounted inside truck cabine)**

**Item number:** 90058578



# We are here for you!

Personal. Capable. Reliable.



Concrete pumps



Mixing plants



Truck mixers



Conveyor belts



Measuring technology



Mixing systems



Concrete recycling systems

## Your powerful partner.

- Delivery and instruction by experienced service technicians
- Worldwide service network with over 90 service centres
- Technical support worldwide
- Fast availability of spare parts
- Professional maintenance and services

Liebherr-Mischtechnik GmbH · Postbox 145 · 88423 Bad Schussenried, Germany  
+49 7583 949-0 · [www.liebherr.com](http://www.liebherr.com) · [info.lmt@liebherr.com](mailto:info.lmt@liebherr.com)

Get your  
contact

