

# ECBN 8872

Integrated fridge-freezer with  
BioFresh and NoFrost

# LIEBHERR

Fully Integrated



Fully Integrated


# ECBN 8872

Integrated fridge-freezer with  
BioFresh and NoFrost

← 91 → 203  
↑  
↓



LIEBHERR



|                                       |   |
|---------------------------------------|---|
| Type                                  | Fixed door / Built-in appliance for integrated use  |
| Total net                             | Fridge: 357 l /<br>Freezer compartment: 114 l   |
| Noise output                          | 39 dB(A)  |
| Climate rating                        | SN-T (+10 °C to +43 °C)   |
| Frequency / voltage                   | 50/60 Hz / 220-240 V~   |
| Connection Rating                     | 2.5 A 189.4 W   |
| Aperture dimensions in cm (h / w / d) | 203.2-207.2 / 91.5 / min. 63.5  |
| Controls                              |   |
| Type of Control                       | MagicEye, touch electronic<br>Digital temperature display for fridge, BioFresh and freezer compartments |
| Door alarm                            | acoustic  |
| Alarm in case of malfunction          | visual and acoustic   |
| Fridge compartment                    |   |
| Forced air system                     | PowerCooling  |
| Door interior                         | Comfort GlassLine storage shelves with bottle holder, egg tray  |
| Bottle shelf                          | Integrated bottle shelf   |
| Interior light                        | LED light column on both sides and LED ceiling lighting<br>LED BioFresh lighting                        |
| BioFresh-Safes                        | 2 BioFresh safes on telescopic rails<br>set variably to either HydroSafe or DrySafe                     |
| Freezer compartment                   |   |
| Star coding                           | ****  |
| Interior light                        | LED lighting  |
| Number of drawers                     | 2 drawers on telescopic rails   |
| IceMaker                              | with fixed 3/4" water connection  |
| Ice cube capacity in 24 hours         | 1.4 kg  |
| Ice cube stock                        | 2.7 kg  |
| Freezing capacity in 24 hours         | 9.5 kg  |
| Temperature rise time                 | 16.9 h  |
| Features                              |   |
| Door hinges                           | Door hinge left fixed/right fixed   |
| Replaceable door seal                 | yes   |
| EAN-Nr.                               | 9005382271254   |

## Features



### BioFresh

Sensitive food also stays fresh for a long time in BioFresh safes. For example, a temperature of close to 0 °C means the Meat & Dairy safe provides the ideal storage conditions for meat or fresh milk. There's also higher humidity in the Fruit & Vegetable safe thanks to the air-tight closure – perfect for fruit and vegetables. The best part is: all optimal adjustments are made in the factory.



### DuoCooling

DuoCooling uses two completely separate refrigeration circuits to ensure no air is exchanged between the fridge and freezer compartments. The food neither dries out nor is odour transferred. This means throwing less away, shopping less often, but saving more and greater enjoyment.




### IceMaker with fixed water connection

Are there any ice cubes left? The IceMaker means you never have to ask. With the IceMaker with fixed water connection, the water comes directly from the tap. The MaxIce function allows you to produce up to 1.5 kg of ice cubes a day. The versatile ice compartment separator lets you choose the amount of ice cubes you want to store; easy and just as you like it.




### NoFrost

When you open your freezer, you want to see frozen food – and certainly not ice and frost. NoFrost protects the freezer compartment from unwanted icing, which consumes a lot of energy and can be expensive. NoFrost means no more tedious and time-consuming defrosting of the freezer compartment, more time for other things – and saving money.



### SuperCool


SuperCool generates greater cold reserves so that freshly stored food can be rapidly chilled. Reversion from the chilling temperature of +2 °C to the normal refrigeration temperature is time or quantity-controlled to help save energy.



### SuperFrost

Want to freeze your fresh food in a way that preserves vitamins? Simply activate SuperFrost – and the freezer compartment will increase its cooling capacity. In this way, the temperature remains constant when you place in unfrozen food. SuperFrost is also energy-saving: if you don't switch off the function manually, it automatically deactivates itself after 65 hours.

## Technical drawing



Products and their details are staged for advertising purposes. Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any damage. (16.08.2024)

