

BCv 1103

CounterTop Cooler




LIEBHERR

BCv 1103

Premium

CounterTop Cooler




home.liebherr.com

BCv 1103	
Total gross/net capacity l (cu.ft.)	3.9/3.2 (110/90)
Exterior dimensions in inches (w/d/h) in mm (w/d/h)	19 5/8"/21 5/8"/32 11/16" 497 / 548 / 830
Interior dimensions in inches (w/d/h) in mm (w/d/h)	16 7/16"/15 7/16"/25 11/16" 417 / 391 / 651
Energy Consumption in 365 days*	510 kWh
Ambient temperature in °C (°F)	50°F to 86°F (10°C to 30°C)
Refrigerant	R600a
Sound rating**	49 dB(A)
Connection rating	2.5 A / 140 W
Frequency/voltage	60 Hz/115V~
Cooling system/Defrost	Forced-air / Manual
Temperature range in °C (°F)	14°F to -13°F (-10 °C to -25 °C)
External cabinet finish	Steel / silver
Door material	Glass door
Internal liner material	Steel
Type of control / Temperature display	Analog control / Internal analogue
Interior Light	LED light column on right
Adjustable shelves	2
Shelf material	Wire shelves plastic -coated
Usable shelf area in inches (mm)	16 7/16"/13 1/8" (417/332)
Shelf load weight in lbs (kg)	66.1 (30)
Handle	Recessed grip
Lock	Recessed
Self-closing door	Yes
Door hinges	Right, reversible on site
Insulation in inches (mm)	1 5/8" (40)
Gross/ net weight in lbs (kg)	95.2/89.7 (43.2/40.7)



Unit Dimensions



Accessories

Scannerrails	7444050
Stacking kit left	9875033
Stacking kit right	9875035

* Measured unloaded at 26.6°C (80°F) ambient temperature with set point of 3.3°C (38°F).

** Noise level – sound power re 1pW – according to IEC 60704-2-14.

LIEBHERR