

EFI 4453

Impulse Sales Chest Freezer




LIEBHERR

EFI 4453

Chest Freezer with Sliding Angle Top

Impulse Sales Chest Freezer




home.liebherr.com



EFI 4453	
Total gross/net capacity l (cu.ft.)	15.0 / 10.7 (426 / 302)
Exterior dimensions in inches (w/d/h) in mm (w/d/h)	57 11/16" / 26 13/16" / 32 1/2" 1465 / 680 / 825
Interior dimensions in inches (w/d/h) in mm (w/d/h)	52 13/16" / 21 1/8" / 19 5/16" 1340 / 535 / 490
Total display area in m ² (TDA)/ft ²	0.79 / 8.50
Energy Consumption in 365 days*	774 kWh
Ambient temperature in °C (°F)	50°F to 86°F (10°C to 30°C)
Refrigerant	R290
Sound rating**	52 dB(A)
Connection rating	2.5 A/245 W
Frequency/voltage	60 Hz/115 - 127V~
Cooling system/Defrost	Static/Manual
Temperature range in °C (°F)	14°F to -11°F (-10°C to -24°C)
Lid material	Sliding glass lid
External cabinet finish	Steel/White
Internal liner material	Steel
Frame color	Grey
Type of control / Temperature display	Analog control/internal analogue
Interior Light	LED
Number of baskets	6
Number of baskets, max. (see picture)	6
Castors	Swivel casters
Lock	can be retrofitted
Insulation in inches (mm)	2 3/8" (60)
Water drain	Yes
Gross/net weight in Kg (lbs)	168.9 / 149.3 (76.6 / 67.7)

Unit Dimensions



Accessories

Basket 10 5/8 in (270 mm)	7111009
Basket 8 1/4 in (210 mm)	7111039
Partition grid 10 5/8 in (270 mm) / High grid	7111031
Partition grid 10 5/8 in (270 mm) / Low grid	7111033
Partition grid 8 1/4 in (210 mm)	7111787
Rubber bumpers, gray	7266707
Cylinder lock for glass sliding lid	7043359
Glass canopy	9809160

* Measured unloaded at 26.6°C (80°F) ambient temperature with set point of 3.3°C (38°F).

** Noise level – sound power re 1pW – according to IEC 60704-2-14.

LIEBHERR

Subject to modification. For up-to-date information, visit home.liebherr.com.