



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Liebherr - LFT50W2HC : LFT50W2HC

Specifications

ENERGY STAR Unique ID:	2323105
ENERGY STAR Partner:	Liebherr-Hausgerate Ochsenhausen GmbH
Brand Name:	Liebherr
Model Name:	LFT50W2HC
Model Number:	LFT50W2HC
Product Type:	High Performance Freezer
Total Volume (cubic feet):	46.79
Lab Grade Refrigerator or Freezer Energy Consumption (kWh/day):	8.59
Peak Temperature Variance (°C):	6.4
Product Form Factor:	Upright
Total Number of Doors:	2
Number of glass doors:	0
Number of solid doors:	2
Door Opening Orientation:	Swinging
Number of Drawers or Shelves:	8
Number of Internal Ports and/or Access Holes:	1
Depth (in):	32.88
Height:	85.06
Width:	56.32
Defrost Type:	Automatic
Refrigerant with GWP:	R-290 (GWP:3 Lower GWP)
Installed Accessories:	None
Date Available on Market:	1514782800000
Date Certified:	1531886400000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Captured On:
06/10/2024