




ELEVATION LEFT ECBN 6156
1 : 10



ELEVATION FRONT ECBN 6156
1 : 10



FLOOR PLAN ECBN 6156
1 : 5



3D View ECBN 6156

ECBN 6156 PremiumPlus BioFresh NoFrost	
Total capacity litre	480
Refrigerator comp. litre	364
BioFresh comp. litre	69
NoFrost Freezer comp. litre	116
Shipping weight kg	177
Climate rating	T
Ice cube output kg / 24 h	1.3
Ice cube stock kg	2.7
Door hinges	right, reversible
Outside dimensions in cm (H / W / D)	203.2-207.2 / 91.5 / at least 62.5
Energy supply Voltage V/Hz	220 – 240 / 50
Electronic control	MagicEye behind-the-door, touch electronic Digital temperature display for fridge and freezer compartments
Permissible Door Weights	Panel dimensions Min. panel thickness = 16 mm Max. panel thickness = 19 mm Maximum panel weight Refrigerator panel = 27 kg Freezer drawer panel = 10 kg

ECBN 6156_136

Scale

As indicated